Setting up TPU for faster processing.

```
%tensorflow_version 2.x
import tensorflow as tf
print("Tensorflow version " + tf. version )
try:
 tpu = tf.distribute.cluster resolver.TPUClusterResolver() # TPU detection
 print('Running on TPU ', tpu.cluster spec().as dict()['worker'])
except ValueError:
 raise BaseException('ERROR: Not connected to a TPU runtime; please see the previous cell in t
tf.config.experimental_connect_to_cluster(tpu)
tf.tpu.experimental.initialize_tpu_system(tpu)
tpu strategy = tf.distribute.TPUStrategy(tpu)
print('Number of replicas:', tpu strategy.num replicas in sync)
    Tensorflow version 2.8.0
    Running on TPU ['10.15.131.250:8470']
    INFO:tensorflow:Deallocate tpu buffers before initializing tpu system.
    INFO:tensorflow:Deallocate tpu buffers before initializing tpu system.
    INFO:tensorflow:Initializing the TPU system: grpc://10.15.131.250:8470
    INFO:tensorflow:Initializing the TPU system: grpc://10.15.131.250:8470
    INFO:tensorflow:Finished initializing TPU system.
    INFO:tensorflow:Finished initializing TPU system.
    INFO:tensorflow:Found TPU system:
    INFO:tensorflow:Found TPU system:
    INFO:tensorflow:*** Num TPU Cores: 8
    INFO:tensorflow:*** Num TPU Cores: 8
    INFO:tensorflow:*** Num TPU Workers: 1
    INFO:tensorflow:*** Num TPU Workers: 1
    INFO:tensorflow:*** Num TPU Cores Per Worker: 8
    INFO:tensorflow:*** Num TPU Cores Per Worker: 8
    INFO:tensorflow: *** Available Device: _DeviceAttributes(/job:localhost/replica:0/task:0/de
    INFO:tensorflow:*** Available Device: DeviceAttributes(/job:localhost/replica:0/task:0/de
    INFO:tensorflow:*** Available Device: _DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: _DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: _DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: _DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: _DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: _DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: _DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: _DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: _DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: _DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: _DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: DeviceAttributes(/job:worker/replica:0/task:0/devic
    INFO:tensorflow:*** Available Device: _DeviceAttributes(/job:worker/replica:0/task:0/devic
```

Introduction ?:




The aim of this project is to create a deep learning algorithum that is able to accurately characterize pneumonia via chest x-rays for pediatric patients. The data set that is being used is coming from Kaggle.com and is being utilized for this project.


→ Process ::

Workflows illustrating the process typically taken for creating a DL project.


Data Preparation




Neural Network Architecture




Fitting the Model



Model Predictions and Evaluation





```
!pip install colorama
!pip install opendatasets
    Collecting colorama
      Downloading colorama-0.4.4-py2.py3-none-any.whl (16 kB)
    Installing collected packages: colorama
    Successfully installed colorama-0.4.4
    Collecting opendatasets
      Downloading opendatasets-0.1.22-py3-none-any.whl (15 kB)
    Requirement already satisfied: click in /usr/local/lib/python3.7/dist-packages (from opence
    Requirement already satisfied: kaggle in /usr/local/lib/python3.7/dist-packages (from oper
    Requirement already satisfied: tqdm in /usr/local/lib/python3.7/dist-packages (from openda
    Requirement already satisfied: certifi in /usr/local/lib/python3.7/dist-packages (from kac
    Requirement already satisfied: urllib3 in /usr/local/lib/python3.7/dist-packages (from kag
    Requirement already satisfied: six>=1.10 in /usr/local/lib/python3.7/dist-packages (from )
    Requirement already satisfied: requests in /usr/local/lib/python3.7/dist-packages (from ka
    Requirement already satisfied: python-slugify in /usr/local/lib/python3.7/dist-packages (f
    Requirement already satisfied: python-dateutil in /usr/local/lib/python3.7/dist-packages (
    Requirement already satisfied: text-unidecode>=1.3 in /usr/local/lib/python3.7/dist-packac
    Requirement already satisfied: idna<3,>=2.5 in /usr/local/lib/python3.7/dist-packages (fro
    Requirement already satisfied: chardet<4,>=3.0.2 in /usr/local/lib/python3.7/dist-packages
    Installing collected packages: opendatasets
    Successfully installed opendatasets-0.1.22
import numpy as np
import pandas as pd
import re
import os
import opendatasets as od
import matplotlib.pyplot as plt
```

from tensorflow import keras

import colorama

from tensorflow.keras.models import Sequential

from tensorflow.keras.layers import Conv2D, MaxPooling2D
from sklearn.model selection import train test split

from tensorflow.keras.preprocessing.image import ImageDataGenerator

from tensorflow.keras.layers import Dense, Dropout, Activation, Flatten

Execution

```
#Importing data set kaggle API for inital data reivew

dataset = 'https://www.kaggle.com/datasets/paultimothymooney/chest-xray-pneumonia'

od.download(dataset)

Please provide your Kaggle credentials to download this dataset. Learn more: <a href="http://bit.ly">http://bit.ly</a>
Your Kaggle username: mdespinoza
Your Kaggle Key: ......

Downloading chest-xray-pneumonia.zip to ./chest-xray-pneumonia
100%| 2.29G/2.29G [00:12<00:00, 200MB/s]
```

▼ Data Review

In this section we will perform an assessment of the dataset without manipulating any of the data.

After importing the data we will need to go through and identify the following:

- · How many data sets are we given?
- · How many classes are there in the dataset?

```
data_dir = './chest-xray-pneumonia'
os.listdir(data_dir)
    ['chest_xray']

# Identifying the number folders
os.listdir('./chest-xray-pneumonia/chest_xray')

    ['val', 'train', 'test', '__MACOSX', 'chest_xray']

# Identifying the categories
os.listdir('./chest-xray-pneumonia/chest_xray/test')
    ['NORMAL', 'PNEUMONIA']


CATEGORIES = ["PNEUMONIA", "NORMAL"]

# Identify the number of images per folder
test_set = ('./chest-xray-pneumonia/chest_xray/test')
train_set = ('./chest-xray-pneumonia/chest_xray/train')
val_set = ('./chest-xray-pneumonia/chest_xray/val')

# function to calculate the number of images in a dataset folder
```


```
def count xray images(title, dataset, color):
    print(color + title)
    pneumonia = len(os.listdir(os.path.join(dataset, 'PNEUMONIA')))
    normal = len(os.listdir(os.path.join(dataset, 'NORMAL')))
    print(f"Pnuemonia = {pneumonia}")
    print(f"Percent Pnuemonia = {pneumonia/(pneumonia + normal)*100:.1f}%")
    print(f"Normal = {normal}")
    print(f"Percent Normal = {normal/(pneumonia + normal)*100:.1f}%")
    print(f"TOTAL = {pneumonia + normal}")
# function to create a bargraph for the number of images per category in any dataset folder
def data visual(title, dataset, y limit):
    pneumonia = len(os.listdir(os.path.join(dataset, 'PNEUMONIA')))
    normal = len(os.listdir(os.path.join(dataset, 'NORMAL')))
    x pos = CATEGORIES
   y = (normal, pneumonia)
   plt.title(title)
   # define color for bars
   c = ("#FF0000", "#ADD8E6")
    # Adding labels, and the title
   plt.xlabel("Category")
    plt.ylabel("Num. of X-Ray Images")
   plt.bar(x pos, y, color = c, ec = "black")
   plt.ylim(y limit)
    for i in range(len(x_pos)):
        plt.text(i, y[i], y[i], ha="center", va="bottom")
    plt.show()
# Calling the function to provide the calculations for each folder and visualization
count xray images("Test Folder", test set, colorama.Fore.BLUE)
data_visual("Test Folder", test_set,[0, 450])
count xray images ("Train Folder", train set, colorama. Fore. GREEN)
data_visual("Train Folder", train_set,[400, 4500])
count_xray_images("Val Folder", val_set, colorama.Fore.RED)
data_visual("Val Folder", val_set,[0, 10])
```

```
Test Folder
     Pnuemonia = 390
     Percent Pnuemonia = 62.5%
     Normal = 234
     Percent Normal = 37.5%
     TOTAL = 624
                              Test Folder
        450
        400
                                             390
        350
     Num. of X-Ray Images
        300
        250
        200
        150
        100
         50
          0
                   PNEUMONIA
                                           NORMAL
                               Category
     Train Folder
     Pnuemonia = 3875
     Percent Pnuemonia = 74.3%
     Normal = 1341
     Percent Normal = 25.7%
     TOTAL = 5216
                               Train Folder
        4500
        4000
                                              3875
        3500
     Num. of X-Ray Images
        3000
        2500
        2000
       1500
        1000
         500
                    PNEUMONIA
                                            NORMAL
                                Category
     Val Folder
     Pnuemonia = 8
     Percent Pnuemonia = 50.0%
# function to display images
def image_xray_batch(path1, path2, title):
    image_list = os.listdir(path1)
    image = path2
    plt.figure(figsize=(5,5))
    for i in range(16):
         plt.subplot(4, 4, i + 1)
         img = plt.imread(os.path.join(image, image list[i]))
         plt.imshow(img, cmap='gray')
         plt.axis('off')
```



```
norm_img = os.listdir("./chest-xray-pneumonia/chest_xray/train/NORMAL")[0]
norm_dir = "./chest-xray-pneumonia/chest_xray/train/NORMAL"
sample_img = plt.imread(os.path.join(norm_dir, norm_img))
plt.imshow(sample_img, cmap='gray')
plt.colorbar()
```

```
plt.title('X Ray Image')
print(f"The max px value = {sample img.max():.1f} and the min = {sample img.min():.1f}")
```



Data Review Summary:

The scope of this project is based on a classification problem that is to asses and predict if a Pediatric patient has pnuemonia or not (2 Categories = Pneumoina, Normal). However in reviewing the training data it is evident that the observations for "Normal" (majority class), is significantly higher than the "Pneumonia" (minority class) observations. Which implies that the dataset is imbalanced and will need to be corrected in the data preparation to avoid any bias towards the "Normal" class.

The process to balanace the data between both classes would be to adjust the weights of each class.

Also, the validation data is very small n=16 and not sufficent for creating a model that we could trust. Thus in the next section we will also split the training and validation data using a standard split of 80:20 to create a larger validation dataset

▼ Data Preparation ¹²/₃₄

```
#accessing the data via googlecloudstorage - which allows for the utliizaiton of Google's TPU

GCS_PATH = 'gs://kds-7d8d6352305c640b4084285616169079853d84b5df55d6be0bc51958'

# paths defined to variables
test_image_folder = (GCS_PATH + '/chest_xray/test/*/*')
train_set_folder = (GCS_PATH + '/chest_xray/train/*/*')
val_set_folder = (GCS_PATH + '/chest_xray/val/*/*')

# Split our trianing dataset into a training and validaiton set via "train_test_split" by a spl
train_files = tf.io.gfile.glob(str(train_set_folder))
train_files.extend(tf.io.gfile.glob(str(val_set_folder)))
```

```
training, validation = train test split(train files, test size=0.2)
# Verification of (n) of train dataset in each category folder created after the split.
normal_n = len([train_files for train_files in training if "NORMAL" in train_files])
print("Normal images: " + str(normal_n))
pneu_n = len([train_files for train_files in training if "PNEUMONIA" in train_files])
print("PNEUMONIA images: " + str(pneu n))
    Normal images: 1063
    PNEUMONIA images: 3122
# creating slices from our defined dataset arrays
# training
data struc train = tf.data.Dataset.from tensor slices(training)
# validation
data_struc_val = tf.data.Dataset.from_tensor_slices(validation)
# test
test_data_struc = tf.data.Dataset.list_files(str(test_image_folder))
# train data
train count = tf.data.experimental.cardinality(data_struc_train).numpy()
print("Training (80%) = " + str(train_count))
# val data
val_count = tf.data.experimental.cardinality(data_struc_val).numpy()
print("Validating (20%) = " + str(val_count))
# test data
test count = tf.data.experimental.cardinality(test data struc).numpy()
print("Test = " + str(test count ))
    Training (80\%) = 4185
    Validating (20\%) = 1047
    Test = 624
# defining data labels as an array
categories_labels = np.array([str(tf.strings.split(item, os.path.sep)[-1].numpy())[2:-1]
                              for item in tf.io.gfile.glob(str(GCS_PATH + '/chest_xray/train/*'
categories labels
    array(['NORMAL', 'PNEUMONIA'], dtype='<U9')</pre>
# Setting dataset to be labeled with categories(0)|(1)
def define label(file path):
    parts = tf.strings.split(file_path, os.path.sep)
    return parts[-2] == "PNEUMONIA"
# defining image size constant
IMAGE\_SIZE = [180, 180]
# Establish unfiormity with the images (3D Tensor) and scaling them down
# to smaller sizes.
```

```
def img conversion(img):
  img = tf.image.decode jpeg(img, channels=3)
  img = tf.image.convert image dtype(img, tf.float32)
 return tf.image.resize(img, IMAGE SIZE)
def path(file path):
    label = define_label(file_path)
    img = tf.io.read_file(file_path)
    img = img conversion(img)
    return img, label
autotune = tf.data.experimental.AUTOTUNE
BATCH SIZE = 16 * tpu strategy.num replicas in sync
train data st = data struc train.map(path, num parallel calls=autotune)
val_data_st = data_struc_train.map(path, num_parallel_calls=autotune)
test_data_st = test_data_struc.map(path, num_parallel_calls=autotune)
test data st = test data st.batch(BATCH SIZE)
# verify data normalization = image shape and label
for image, label in train data st.take(1):
 print("Image shape: ", image.numpy().shape)
 print("Label: ", label.numpy())
    Image shape: (180, 180, 3)
    Label: True
def prep_data_pre_training(ds, cache=True, shuffle_buffer_size=1000):
    # This is a small dataset, only load it once, and keep it in memory.
    # use `.cache(filename)` to cache preprocessing work for datasets that don't
    # fit in memory.
    if cache:
        if isinstance(cache, str):
            ds = ds.cache(cache)
        else:
            ds = ds.cache()
    ds = ds.shuffle(buffer_size=shuffle_buffer_size)
    # Repeat forever
   ds = ds.repeat()
   ds = ds.batch(BATCH SIZE)
    # `prefetch` lets the dataset fetch batches in the
    # background while the model is training.
    ds = ds.prefetch(buffer_size=autotune)
```

```
train_data_st = prep_data_pre_training(train_data_st)
val_data_st = prep_data_pre_training(val_data_st)

image_batch, label_batch = next(iter(train_data_st))

# Adjusting the weights for both classes to compensate for the
# data not being equal
weight_for_0 = (1 / normal_n)*(train_count)/2.0
weight_for_1 = (1 / pneu_n)*(train_count)/2.0

class_weight = {0: weight_for_0, 1: weight_for_1}

print('Normal Weight: {:.2f}'.format(weight_for_0))
print('Pneumonia Weight: {:.2f}'.format(weight_for_1))

Normal Weight: 1.97
Pneumonia Weight: 0.67
```

▼ Neural Network Architecture ¾

```
def build model():
    # The model is going to consist of a sequence of layers as
   # noted by tf.keras.Squential
   model = tf.keras.Sequential([
        tf.keras.Input(shape=(IMAGE_SIZE[0], IMAGE_SIZE[1], 3)), # expected input shape
        # Conv2D = convolution layer
        # setting filters (16) and kernel size (3)
        # padding defined as 'same' which implies the preservation
        # of spatial dminensions
        # with activation Rectified Linear Unit = 0 or greater is the ouput
        tf.keras.layers.Conv2D(16, 3, activation='relu', padding='same'),
        tf.keras.layers.Conv2D(16, 3, activation='relu', padding='same'),
        # Maxpool2D - pooling of max values
        tf.keras.layers.MaxPool2D(),
        tf.keras.Sequential([
        tf.keras.layers.SeparableConv2D(32, 3, activation='relu', padding='same'),
        tf.keras.layers.SeparableConv2D(32, 3, activation='relu', padding='same'),
        tf.keras.layers.BatchNormalization(),
        tf.keras.layers.MaxPool2D()]),
        tf.keras.Sequential([
        tf.keras.layers.SeparableConv2D(64, 3, activation='relu', padding='same'),
        tf.keras.layers.SeparableConv2D(64, 3, activation='relu', padding='same'),
        tf.keras.layers.BatchNormalization(),
        tf.keras.layers.MaxPool2D()]),
        tf.keras.Sequential([
        tf.keras.layers.SeparableConv2D(128, 3, activation='relu', padding='same'),
        tf.keras.layers.SeparableConv2D(128, 3, activation='relu', padding='same'),
```

```
tf.keras.layers.MaxPool2D()]),
        tf.keras.layers.Dropout(0.2),
        tf.keras.Sequential([
        tf.keras.layers.SeparableConv2D(256, 3, activation='relu', padding='same'),
        tf.keras.layers.SeparableConv2D(256, 3, activation='relu', padding='same'),
        tf.keras.layers.BatchNormalization(),
        tf.keras.layers.MaxPool2D()]),
        # prevents over fitting the model
        tf.keras.layers.Dropout(0.2),
        # reshape the input
        tf.keras.layers.Flatten(),
        tf.keras.Sequential([
        # Dense implies adding a layer of neurons
        tf.keras.layers.Dense(512, activation='relu'),
        # normalizing the data
        tf.keras.layers.BatchNormalization(),
        tf.keras.layers.Dropout(0.7)]),
        tf.keras.Sequential([
        tf.keras.layers.Dense(128, activation='relu'),
        tf.keras.layers.BatchNormalization(),
        tf.keras.layers.Dropout(0.5)]),
        tf.keras.Sequential([
        tf.keras.layers.Dense(64, activation='relu'),
        tf.keras.layers.BatchNormalization(),
        tf.keras.layers.Dropout(0.3)]),
        tf.keras.layers.Dense(1, activation='sigmoid')
    1)
    return model
# Distribute training across the TPU
with tpu strategy.scope():
   model = build model()
    # defining metrics to be as outputs in the neural network
   METRICS = ['accuracy', # fration of predictions the model gets correct
        # measure what proportion of "+" id's were correct
        tf.keras.metrics.Precision(name='precision'),
        # measure what prportion of actual positives weere id correctly
        tf.keras.metrics.Recall(name='recall')]
    #setting up compiler
    model.compile(
        # Applying gradient descent algorithm
        optimizer='adam',
        # cross entropy loss for binary (0 or 1 )
        # classifications = NORMAL or PNEUMONA
```

tf.keras.layers.BatchNormalization(),

→ Fitting the Model X

```
EPOCHS = 15
#initial - 1st model fit training
history = model.fit(
 train data st,
 steps per epoch= train count // BATCH SIZE, # num of sample / batch size
 epochs=EPOCHS,
 validation data=val data st,
 validation steps= val count // BATCH SIZE, # num of sample / batch size
 class weight=class weight, # data imbalance correction by weights
)
 Epoch 1/15
 Epoch 2/15
 Epoch 3/15
 Epoch 4/15
 Epoch 5/15
 Epoch 6/15
 Epoch 7/15
 Epoch 8/15
 Epoch 9/15
 Epoch 10/15
 Epoch 11/15
 Epoch 12/15
 Epoch 13/15
 Epoch 14/15
 Epoch 15/15
 model.summary()
 Model: "sequential_7"
```

```
max pooling2d (MaxPooling2D (None, 90, 90, 16)
     sequential (Sequential) (None, 45, 45, 32)
                                                        2160
     sequential_1 (Sequential) (None, 22, 22, 64)
                                                        7392
     sequential 2 (Sequential) (None, 11, 11, 128)
                                                        27072
     dropout (Dropout)
                               (None, 11, 11, 128)
     sequential 3 (Sequential) (None, 5, 5, 256) 103296
                              (None, 5, 5, 256)
     dropout_1 (Dropout)
     flatten (Flatten)
                              (None, 6400)
     sequential_4 (Sequential) (None, 512)
                                                        3279360
     sequential 5 (Sequential) (None, 128)
                                                        66176
     sequential_6 (Sequential) (None, 64)
                                                        8512
     dense 3 (Dense)
                               (None, 1)
                                                        65
    ______
    Total params: 3,496,801
    Trainable params: 3,494,433
    Non-trainable params: 2,368
# Save keras model so that model may be used later to continue training the data
# from the saved state.
checkp_cb = tf.keras.callbacks.ModelCheckpoint("pneumonia_model.h5",
                                            save_best_only=True)
# training will terminate if loss does not improve after "5" epochs
early stop cb = tf.keras.callbacks.EarlyStopping(patience=5,
                                              restore best weights=True)
# defining learning rate decay to be applied accross the 2nd model fit training
# until local minima is obtained
def exponential decay(lr0, s):
   def exponential_decay_fn(epoch):
       return lr0 * 0.1 **(epoch / s)
   return exponential_decay_fn
exponential_decay_fn = exponential_decay(0.01, 20)
lr_scheduler = tf.keras.callbacks.LearningRateScheduler(exponential_decay_fn)
#second - 2nd model fit training
history = model.fit(
   train data st,
   steps per epoch= train count // BATCH SIZE,
   # increased batch size as callbacks has been added to stop training after
   # defined criteria has been achieved.
   epochs= 100,
```

```
validation_data=val_data_st,
validation_steps= val_count // BATCH_SIZE,
class_weight=class_weight,
callbacks=[checkp_cb, early_stop_cb, lr_scheduler]
```

```
Epoch 1/100
Epoch 2/100
Epoch 3/100
Epoch 4/100
Epoch 5/100
Epoch 6/100
Epoch 7/100
Epoch 8/100
Epoch 9/100
Epoch 10/100
Epoch 11/100
Epoch 12/100
Epoch 13/100
Epoch 14/100
Epoch 15/100
Epoch 16/100
Epoch 17/100
Epoch 18/100
Epoch 19/100
Epoch 20/100
Epoch 21/100
Epoch 22/100
Epoch 23/100
Epoch 24/100
Epoch 25/100
Epoch 26/100
Epoch 27/100
Epoch 28/100
32/32 [========================] - 48s 2s/step - loss: 0.0293 - accuracy: 0.9885 - r
```


model.summary()

Model: "sequential_7"

Layer (type)	Output Shape	Param #
conv2d (Conv2D)	(None, 180, 180, 16)	
conv2d_1 (Conv2D)	(None, 180, 180, 16)	2320
<pre>max_pooling2d (MaxPooling2D)</pre>	(None, 90, 90, 16)	0
sequential (Sequential)	(None, 45, 45, 32)	2160
sequential_1 (Sequential)	(None, 22, 22, 64)	7392
sequential_2 (Sequential)	(None, 11, 11, 128)	27072
dropout (Dropout)	(None, 11, 11, 128)	0
sequential_3 (Sequential)	(None, 5, 5, 256)	103296
dropout_1 (Dropout)	(None, 5, 5, 256)	0
flatten (Flatten)	(None, 6400)	0
sequential_4 (Sequential)	(None, 512)	3279360
sequential_5 (Sequential)	(None, 128)	66176
sequential_6 (Sequential)	(None, 64)	8512
dense_3 (Dense)	(None, 1)	65
Total params: 3,496,801 Trainable params: 3,494,433 Non-trainable params: 2,368		=======

```
# plotting training results
fig, ax = plt.subplots(1, 3, figsize=(20, 5))
ax = ax.ravel()

for i, met in enumerate(['precision', 'recall', 'accuracy',]):
    ax[i].plot(history.history[met])
    ax[i].plot(history.history['val_' + met])
    ax[i].set_title('Model {}'.format(met))
    ax[i].set_xlabel('epochs')
    ax[i].set_ylabel(met)
    ax[i].legend(['train', 'val'])
```



Model Predictions

```
loss, acc, prec, rec = model.evaluate(test data st)
```

Summary /

- 1. The model overall has an 98.95% accuracy and a loss of 3.55%.
 - model may be improved by adjusting the CNN (Fine-tunning)
- 2. The model predicition accuracy is 80.29% which is less than both the training and validation accuracy. (Can be improved!)
- 3. The model prediction Recall is larger than the precision. Precision may be imrpoved by adding more data that is of high quality.
 - 76.33 % of total class outcomes was predicted correctly
- 4. To improve the results there are two options:
 - Fine-tunning of the CNN and re-training it.
 - · Capture more data

Sources



Acknowledgements Data: https://data.mendeley.com/datasets/rscbjbr9sj/2

License: CC BY 4.0

Citation: http://www.cell.com/cell/fulltext/S0092-8674(18)30154-5

Kaggle Project: https://www.kaggle.com/datasets/paultimothymooney/chest-xray-pneumonia

Other Kaggle opensource review on dataset:

- 1. https://www.kaggle.com/code/aakashnain/beating-everything-with-depthwise-convolution
- 2. https://www.kaggle.com/code/amyjang/tensorflow-pneumonia-classification-on-x-rays
- 3. https://towardsdatascience.com/deep-learning-for-detecting-pneumonia-from-x-ray-images-fc9a3d9fdba8
- 4. https://towardsdatascience.com/detecting-covid-19-induced-pneumonia-from-chest-x-rays-with-transfer-learning-an-implementation-311484e6afc1

Colab TPU Use With Flowes modeling:

https://colab.research.google.com/notebooks/tpu.ipynb#scrollTo=LtAVr-4CP1rp

Webiste:

 $\underline{https://colab.research.google.com/notebooks/snippets/accessing_files.ipynb\#scrollTo=z1_FuDjAozF1}$

X